



DAVID
DAVID MEDICAL

>> JH20-1 Series

Transcutaneous Jaundice Detector

JH20-1A



JH20-1B



JH20-1C



NINGBO DAVID MEDICAL DEVICE CO., LTD

Endless Love keeps the new life safe

**JH20-1 Series
Transcutaneous
Jaundice Detector**



Technical specification:

Model	JH20-1A	JH20-1B	JH20-1C
Checking methods	Compare the blue and green light	Compare the blue and green light	Compare the blue and green light
Indication	Two-figure LED digital display	Three-figure LED digital display(one decimal digits)	Three-figure LCD digital display(one decimal digits)
Readable	Switch the measured value into the bilirubin value according to the comparative table	Test the bilirubi value directly by the forehead of baby	Test the bilirubi value directly by the forehead of baby
Average counting	No	No	Can calculate the measured average values for 2-5times automatically
Test deviation	00~20±1 21~40±2	0~15±1 16~25±1.5	0~15±1 16~25±1.5
Light resource	Xenon-flashing light	Xenon-flashing light	Xenon-flashing light
Power supply	Nickel xenon battery No.5 of 4 pieces	Nickel xenon battery group 4.8V	Nickel xenon battery group 4.8V
Prepare time of start-up	<10s	<10s	<5s
Charger	Input 220V 50Hz Output 6V 0.3ADC	Input 220V 50Hz Output 6V 0.3ADC	Input 220V 50Hz Output 6V 0.3ADC
Checkout	White screen indicates 00 or 01 Yellow screen indicates 20.0±1	White screen indicates 00 or 01 Yellow screen indicates 20.0±1	White screen indicates 00 or 01 Yellow screen indicates 20.0±1
Weight	About 340g	About 245g	About 190g
Volume	171mm x 74mm x 41mm	163mm x 66mm x 37mm	161mm x 53mm x 30mm